




PODO LINE

TECN**DENT**
INNOVATIVE MEDICAL SOLUTIONS





Tecnodent Podiatry chairs line is designed to satisfy every health care provider's need in the fields of podiatry, cosmetics and general practice.

Our Podiatry line products provide the right blend of ergonomics, highest functionality and a reassuring and fine design, round shaped lines, anatomical seamless upholsteries for perfect cleaning and hygiene.

Linda3 PODO is a podiatric chair equipped with electronically controlled movement, independently from height, backrest or leg rest inclination. It offers programmable treatment positions among which Trendelenburg.

Two independent leg rests are included with the **Linda3 PODObis**, which allow the doctor to choose the most proper positioning of the patient's foot and give the professional more freedom and place to operate on one single foot.

Each leg rest can be easily regulated in length, angle and tilting.

Tecnodent **PODOLight** chair's peculiarity is the one-piece backrest, without headrest and fully articulating armrests, for a better accessibility for the patient.

Two independent leg rests provide a functional treatment to the patient.

Linda3 PODO

Multifunction **Linda3 PODO** is equipped with:

- Electromechanical up/down movements, tilting backrest, adjustable leg rest;
 - Anatomical backrest, with lateral arm support wings, to enhance the patient's comfort;
 - Ergonomic PODO headrest;
 - Fully left and right symmetry, for both left and right handed professionals;
 - Mobile high precision pedal control;
 - Emergency chair stop onto both backrest and leg rest;
 - Automatic return to "0" position after the treatment;
 - 3 programmable operating positions;
 - Minimum height 50 cm, which allows easy access also for elder and disabled patients;
- Anatomical upholstery without seams.

Options Linda3 PODO:

- Inox basin with support;
- Magnetic leg rest pillow.

Linda3 PODObis

Multifunction **Linda3 PODObis** is equipped with.

- Electromechanical up/down movements, tilting backrest;
 - Anatomical backrest, with lateral arm support wings, to enhance the patient's comfort;
 - Ergonomic PODO headrest;
 - Fully left and right symmetry, for both left and right handed professionals;
 - Mobile high precision pedal control;
 - Emergency chair stop onto both backrest and leg rest;
 - Independent and manually length, angle and inclination adjustable leg-rests;
 - Anatomical upholstery without seams;
 - Automatic return to "0" position after the treatment;
- Automatic memory of last operating position;
 - Minimum height 50 cm, which allows easy access also for elder and disabled patients.

Option Linda3 PODObis:

- Electronic backrest/leg rest synchronism to reach Trendelenburg position;
- 3 programmable operating positions.



PODOLight

PODOLight is a lean and compact chair, suited for small spaces and for those who want the perfect combination of easy and safe functionality with a competitive and aggressive price.

Two fully articulating armrests help to easily access to the chair. Each leg section can be manually regulated in length, in widening and in inclination. Pedal foot control provides precise adjustability of the chair. Trendelenburg position is reached through electronic backrest/leg rest synchronism. Tiny wheels allow small and precise movements of the chair.





**PODO headrest
(for Linda3 PODO and PODObis)**

Anatomical headrest for a perfect cervical support.



UNI headrest (OPTION, for Linda3 PODO and PODObis)

Multiarticulating headrest for the patient's perfect comfort and positioning.

Instrument tray (OPTION, for Linda3 PODO and PODObis)

Functional instrument tray at the side of the seat.



Right and/or left armrest (OPTION, for Linda3 PODO and PODObis)

Rotating armrests help the patient to get on/off the chair.

Infrared control (OPTION, for Linda3 PODObis)

Remote control which allows to choose and set the chair position.



Mobile foot control "Easy Touch"

Easy and functional foot control which provides different precise treatment positions and doctor's hands-free control.

Magnetic leg rest pillow (OPTION, for Linda3 PODO)

Ergonomic pillow for a comfortable treatment.



Stainless steel basin (OPTION, for Linda3 PODO)

Stainless steel basin that can be easily removed as necessary.

Standard base

Compact polyurethane base which guarantees stability and free room for professional's movements.



Rotating base (OPTION, for Linda3 PODO and PODObis)

It allows a 325° rotation around the chair's longitudinal axis.

TECHNICAL FEATURES		Linda3 PODO		Linda3 PODObis		PODolight
		without rotation	with rotation	without rotation	with rotation	
Maximum height with extended leg rest	mm	1040	1040	1040	1040	1040
Minimum height	mm	500	500	500	500	608
Width	mm	660	610	660	610	655
Rotation angle	°	-	325	-	325	-
Length	mm	1900	1900	1900	1900	1848
Gross weight	kg	145	145	145	145	90
Net weight	kg	115	115	115	115	74
Maximum lifting weight	kg	210	210	210	210	135
Longitudinal sliding of the leg rests		I	I	I	I	I
Lifting angle of the leg rests		B	B	B	B	B
Divarication angle of the legrests	mm			150	150	150
Insulation class	°			100	100	100
Parts type	°			30	30	30